

medw|ork®
for mindful people

CrossSnare ZERO

Cold polyp resection in a new dimension

www.medwork.com

Cold polyp resection
in a new dimension

CrossSnare ZERO

The new dedicated cold snare

Perfect resection thanks to unique tube design

Multiple snare diameters for diminutive polyps

Optimal for the implementation of the ESGE guideline*

#exact

NEW

Great effect on small polyps

Less effort and perfect cutting performance - that's what CrossSnare ZERO offers you thanks to its unique reinforced spiral tube design.

Equipped with extra thin cutting wire and available in two different sizes, CrossSnare ZERO is the perfect choice for polyps up to 9 mm and supports you optimally in implementing the ESGE guideline*.

Product Code	Order No.	Characteristics	ø Snare	Length	ø Working Channel	Unit
② POL1-Z2-07-23-220	502569	CWT, Kite	7 mm	220 cm	2.8 mm	10
② POL1-Z2-12-23-220	502570	CWT, Kite	12 mm	220 cm	2.8 mm	10

* Ferlitsch Monika et al.: „Colorectal polypectomy and endoscopic mucosal resection (EMR): European Society of Gastrointestinal Endoscopy (ESGE) Clinical Guideline“, Endoscopy 2017.